



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 09/98/MW	Status 09/98/MW															
Substance code:		Page : 1 (3)															
1. Identification of the substance / preparation and company Product details Trade Name: Arctic White (Tooth Whitening Gel) Supplier details : PSP Dental Company Limited 3-5 Dylan Road Belvedere Kent DA17 5QS U.K Information provided by: PSP Dental Company Limited Emergency telephone number: 020-8311-7337 Emergency fax number: 020-8310-0920																	
2. Composition/information on ingredients Chemical characterization Hazardous Components <table border="1"><thead><tr><th>Ingredient</th><th>%</th><th>CAS No.</th><th>S Phrase</th><th>R Phrase</th></tr></thead><tbody><tr><td>Organic Peroxide</td><td>20 - 40</td><td>124-43-6</td><td>S24, 25</td><td></td></tr><tr><td>Triethanol Amine</td><td>4 - 8</td><td>102-71-6</td><td>S24, 25</td><td></td></tr></tbody></table>			Ingredient	%	CAS No.	S Phrase	R Phrase	Organic Peroxide	20 - 40	124-43-6	S24, 25		Triethanol Amine	4 - 8	102-71-6	S24, 25	
Ingredient	%	CAS No.	S Phrase	R Phrase													
Organic Peroxide	20 - 40	124-43-6	S24, 25														
Triethanol Amine	4 - 8	102-71-6	S24, 25														
3. Hazards identification Unlikely to cause harmful effects under normal conditions of handling and use. Practice good hygiene. Wear gloves, overalls.																	
4. First aid measures After contact with skin: Wash skin first with water and then, soap and water. If in doubt seek medical help. After contact with eyes: Irritant, wash open eye with of clean water for 15 minutes. Seek medical help. After ingestion: Seek medical advice. After inhalation: N/A																	



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 01/98/HLS	Status 01/98/HLS
Substance code:		Page : 2 (3)
5. Fire-fighting measures: May evolve toxic fumes in fire. Suitable extinguishing media: foam, dry powder, carbon dioxide		
6. Accidental release measures: Environmental precautions Prevent entry into drains and water courses. Dispose as indicated by Local Authority Regulations. Methods for cleaning up / taking up: Wash affected skin areas with soap and water.		
7. Handling and storage Handling: Good General Ventilation. No special precautions are necessary. Follow instructions enclosed in package and use only for purposes directed Avoid contact with eyes and skin. Avoid creating dusts. Storage: Store in a cool dry place. Avoid direct heat and sunlight.		
8. Exposure controls / personal protection Hygiene measures Remove soiled or soaked clothing immediately. Wash hands before breaks and after work. Respiratory protection: N/A Hand protection: rubber gloves. Eye protection: Goggles. Body protection: Overalls and boots		
9. Physical and chemical properties Appearance Form : Gel Colour : Transparent Odour : Vanilla Data relevant to safety Flash point : N/A Melting Point Range: N/A Ignition temperature N/A Density : N/A Explosive Limits : N/A Solubility in water: N/A Boiling Point / Range : N/A Vapour Pressure: N/A		



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 01/98/HLS	Status 01/98/HLS
Substance code:		Page : 3 (3)
10. Stability and reactivity		
Thermal decomposition / Conditions to avoid Stable in normal conditions		
Hazardous reactions / Substances to avoid : Avoid Heating		
11. Toxicological information		
There is no Toxicological data available for this product. There is no known health risk if used properly. Prolonged skin and eye contact should be avoided. Remarks: None		
12. Ecological information		
Further ecological information Remarks: Believed not to be dangerous to the environment if disposed with care and attention.		
13. Disposal considerations		
Product Substance via authorized waste disposal contractor to an approved waste disposal site, observing all local and national regulations.		
14. Transport information		
UN No: N/A	IMDG Class: N/A	
IMO: N/A	Packaging Group: N/A	
IATA: N/A	Packaging Group: N/A	
Correct Technical Name: N/A ADR/RID: N/A Not regulated for transport: N/A		
15. Regulatory information		
Hazard symbols:		
R phrases:		
S phrases:	S25: Avoid contact with eyes S24: Avoid contact with skin	
Other Regulations: N/A		
16. Other information		
This Information is given in good faith, being based on the latest knowledge available to P.S.P.Dental Co. Ltd. No known relevant information has been omitted from the product information, risk information and safety precaution sheets. The user should not assume therefore, the product is suitable for any abnormal use. If the information provided is insufficient to ensure safety in any particular application, please contact P.S.P.Dental Co. Ltd prior to use for further advice.		