



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 09/98/MW	Status 09/98/MW
Substance code:		Page : 1 (3)
1. Identification of the substance / preparation and company		
Product details		
Trade Name: Ceram II VLC Powder		
Supplier details :		
PSP Dental Company Limited		Information provided by: PSP Dental Company Limited
3-5 Dylan Road		
Belvedere		
Kent DA17 5QS		
Emergency telephone number: 020-8311-7337		
Emergency fax number: 020 8310 0920		
2. Composition/information on ingredients		
Chemical characterization		
Hazardous Components		
Ingredient	%	CAS No. S Phrase R Phrase
Barium, Strontium alumino silicate fillers	< 85.00	
Dried polyacrylic acid	< 6.00	
Miscellaneous inorganic fillers	< 15.00	
Benzoyl peroxide	< 0.5	
Non Hazardous Components	%	CAS No.
Ingredient		
3. Hazards identification		
OSHA Permissible Exposure Limit	N/A	
Carcinogen - NIP Programme	N/A	
ACGIH threshold limit value	N/A	
4. First aid measures		
After contact with skin:	Rinse with plenty of water, seek medical help if discomfort continues.	
After contact with eyes:	Rinse with plenty of water. Seek medical help if discomfort continues.	
After ingestion:	Seek medical help immediately	
After inhalation:	Seek medical help immediately	



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 01/98/HLS	Status 01/98/HLS
Substance code:		Page : 2 (3)
5. Fire-fighting measures: Suitable extinguish media Use extinguishing method and media suitable to surrounding fire. Water, Water Spray, Dry Powder Foam CO ²		
6. Accidental release measures: Environmental precautions N/A Methods for cleaning up / taking up: Spill Response: Sweep, then mop with water, do not use bare hands. Watch out for dust hazards.		
7. Handling and storage Handling: Handle with gloves Storage: The container should be stored in normal dry place with lid firmly closed, since material is moisture sensitive.		
8. Exposure controls / personal protection Hygiene measures Remove soiled or soaked clothing immediately. Wash hands before breaks and after work. Respiratory protection: N/A Hand protection: rubber gloves. Eye protection: N/A Body protection: None		
9. Physical and chemical properties Appearance Form : Fine Powder Colour : Varying Colour Odour : None Data relevant to safety Flash point : N/A Melting Point Range: N/A Ignition temperature N/A Density : N/A Explosive Limits : N/A Solubility in water: 10-WT% (Dry Acid Content) Boiling Point / Range : N/A Vapour Pressure:		



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 01/98/HLS	Status 01/98/HLS
Substance code:		Page : 3 (3)
10. Stability and reactivity Thermal decomposition : None Conditions to avoid : >25EC (77EF) and Reducing Agents, e.g: Amines Hazardous reactions : Eugenol Products Substances to avoid: Water		
11. Toxicological information Remarks: LD 50 7710 mg/kg of Benzoyl Peroxide LD 50 > 5000 mg/kg of Barium, strontium aluminosilicate fillers LD 50 > 5000 mg/kg of Dried poly (acrylic acid) LD 50 > 10000 mg/kg of Miscellaneous inorganic fillers		
12. Ecological information Further ecological information Remarks		
13. Disposal considerations Product : Do not flush in to surface water or sanitary sewer system		
14. Transport information UN No: IMDG Class: IMO: Packaging Group: IATA: Packaging Group: Correct Technical Name: ADR/RID: Not regulated for transport		
15. Regulatory information Hazard symbols: R phrases: S phrases: S25: Avoid contact with eyes Other Regulations:		
16. Other information This Information is given in good faith, being based on the latest knowledge available to P.S.P.Dental Co. Ltd. No known relevant information has been omitted from the product information, risk information and safety precaution sheets. The user should not assume therefore, the product is suitable for any abnormal use. If the information provided is insufficient to ensure safety in any particular application, please contact P.S.P.Dental Co. Ltd prior to use for further advice.		