



Dental Company Limited

3/5 DYLAN ROAD  
BELEVEDERE  
KENT  
DA17 5QS

**HEALTH AND SAFETY DATA SHEET**

Product no. Refit Soft/Hard Powder	Version 09/98/MW	Status 04/04/MO
Substance code: 10552		Page : 1 (3)
<b>1. Identification of the substance / preparation and company</b>		
Product details		
Trade Name: Refit Soft / Hard Powder		
Supplier details :		
PSP Dental Company Limited		Information provided by: PSP Dental Company Limited
3-5 Dylan Road		
Belvedere		
Kent DA17 5QS		
U.K		<b>Emergency telephone number: 020-8311-7337</b>
		<b>Emergency fax number: 020 8310 0920</b>
<b>2. Composition/information on ingredients</b>		
Chemical characterization		
Hazardous Components		
Ingredient	%	CAS No. S Phrase R Phrase
Ethyl Methacrylate Polymers	95.00	97-63-2
CC Polymer	4.00	S22,37 R43
Pigments	1.00	
<b>3. Hazards identification</b>		
<b>Unlikely to cause harmful effects under normal conditions of handling and use.</b>		
<b>Practice good hygiene. Wear gloves, overalls and mask.</b>		
<b>4. First aid measures</b>		
After contact with skin:		
Mild irritant, prolonged contact may cause defatting and dermatitis. Wash skin first with water and then, soap and water. If in doubt seek medical help.		
After contact with eyes:		
Mild irritant, wash open eye with of clean water for 10 minutes. Seek medical help.		
After ingestion: Ensure clear airways. Rinse out mouth, in case of ingestion seek medical attention.		
Seek medical advice.		
After inhalation:		
Remove to fresh air. If breathing is irregular or stopped, give artificial respiration. Give nothing by mouth. Seek medical help.		



