



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no. Varnish	Version 01/98/HLS	Status 01/98/HLS
Substance code:		Page : 1 (3)
1. Identification of the substance / preparation and company		
Product details		
Trade Name: Varnish		
Supplier details :		
PSP Dental Company Limited 3-5 Dylan Road Belvedere Kent DA17 5QS		Information provided by: PSP Dental Company Limited
		Emergency telephone number: 0208-311-7337 Emergency fax number: 0208-310 -0920
2. Composition/information on ingredients		
Chemical characterization		
Hazardous Components		
Ingredient	%	CAS No. S Phrase R
Phrase		
Acetone	44.0	67-64-1
Ether	28.0	60-29-7
Cyclohexane	11	110-80-7
Cellulose	3	
Non Hazardous Components		
Ingredient	%	CAS No.
Natural Gum	1 2	
3. Hazards identification		
Unlikely to cause harmful effects under normal conditions of handling and use. Practice good hygiene. Wear gloves, overalls.		
4. First aid measures		
after contact with skin Wash with copious amounts of running water and soap. Seek medical attention if irritation persists.		
after contact with eyes Wash with copious amounts of running water . If discomfort persists seek medical attention without delay.		
after ingestion Rinse mouth thoroughly with water. Give plenty of water to drink, allow the patient to rest and keep warm. Seek medical advice.		
after inhalation Remove subject from the area , allow subject to rest and keep warm.		



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no. Varnish	Version 01/98/HLS	Status 01/98/HLS
Substance code:		Page : 2 (3)
5. Fire-fighting measures Suitable extinguishing media Carbon Dioxide, dry powder, alcohol resistant foam sand, earth or water fog. Do not use water in a jet.		
6. Accidental release measures Environmental precautions Shut off sources of ignition , Wear appropriate protective clothing. Absorb on an inert material, transfer To containers for disposal according to regulations. Wash site with water and detergent. Methods for cleaning up / taking up As above. Equipment can be wiped clean with a Varnish solvent.		
7. Handling and storage Handling Keep away from all sources of ignition. No smoking. Beware of sparks caused by friction or static electricity. Use in a well ventilated place. Wear safety glasses and impervious gloves. Storage The product should not be stored in hot areas. Keep away from direct light.		
8. Exposure controls / personal protection Hygiene measures Remove soiled or soaked clothing immediately. Wash hands before breaks and after work. Respiratory protection :Type approved PPE for dust if required Hand protection :Gloves Eye protection :Goggles Body protection :Overalls and boots		
9. Physical and chemical properties Appearance Form :Liquid Colour :Clear/Pale Yellow Odour :Resinous Odour Data relevant to safety Flash point 60Ec Ignition temperature Explosive Limits Boiling Point / Range Melting Point Range Density Solubility in water Vapour Pressure		



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no. Varnish	Version 01/98/HLS	Status 01/98/HLS										
Substance code:		Page : 3 (3)										
10. Stability and reactivity Thermal decomposition / Conditions to avoid Stable. Avoid all sources of ignition. Hazardous reactions / Substances to avoid Avoid oxidising agents, reducing agents.												
11. Toxicological information There is no Toxicological data available for this product. There is no known health risk if used properly. Prolonged skin and eye contact should be avoided. Remarks												
12. Ecological information Further ecological information Remarks No environmental hazard is anticipated provided that the material is handled and disposed of with due care and attention.												
13. Disposal considerations Product Container and substance via authorised waste disposal contractor to an approved waste disposal site observing All local and National regulations.												
14. Transport information <table border="0" style="width: 100%;"> <tr> <td style="width: 50%;">UN No: 1263</td> <td style="width: 50%;">IMDG Class:</td> </tr> <tr> <td>IMO:</td> <td>Packaging Group:</td> </tr> <tr> <td>IATA:</td> <td>Packaging Group:</td> </tr> <tr> <td>Correct Technical Name:</td> <td></td> </tr> <tr> <td>ADR/RID:</td> <td></td> </tr> </table>			UN No: 1263	IMDG Class:	IMO:	Packaging Group:	IATA:	Packaging Group:	Correct Technical Name:		ADR/RID:	
UN No: 1263	IMDG Class:											
IMO:	Packaging Group:											
IATA:	Packaging Group:											
Correct Technical Name:												
ADR/RID:												
15. Regulatory information Hazard symbols R phrases R 11 Highly Flammable S phrases S9-16-33-23 Keep container in a well ventilated place. Keep away from sources of ignition No smoking. Take precautionary measures against static discharge. Do not breathe vapour. Other regulations												
16. Other information This Information is given in good faith, being based on the latest knowledge available to P.S.P.Dental Co. Ltd. No known relevant information has been omitted from the product information, risk information and safety precaution sheets. The user should not assume therefore, the product is suitable for any abnormal use. If the information provided is insufficient to ensure safety in any particular application, please contact P.S.P.Dental Co. Ltd prior to use for further advice.												